


- **Compatible with IPV devices**  
(except IPV-HC)

- **Features**

- "CE" certified product; FDA registered
- Validated cleaning and disinfection process
- Safety valve designed to protect against risk of internal overpressure
- Harness connections are specially designed to reduce risk and facilitate use
- Instructions for Use & Quick Reference Guide included



- **Product**

Disposable, single patient use  
Non-sterile, ready-to-use from package  
Adult, pediatric, and neonate  
Duration of use: up to six (6) months

- **Specifications**

Overpressure safety valve opens at	50cm H2O internal pressure
Operates at frequencies of	60 – 600 cycles/min
Accommodates a Flowrate up to	170L/min
ISO 5356 Compliant connection port(s)	15mm ID & 22mm OD

- **Materials**

Medical grade plastic; compliant with biocompatibility requirements  
Latex free

- **Standard Packaging**

Units / Bag / Box	1 Unit / Bag / Box
Unit / Case	10 Units / Case

*For further product information, please visit  
CircuitLife® website at [www.circuit-life.com](http://www.circuit-life.com)*

