



AEROCHAMBER PLUS® Z STAT® aVHC

The innovative AEROCHAMBER® brand of valved holding chamber (aVHC) has been trusted by hospitals for over 35 years.

The AEROCHAMBER PLUS® Z STAT® aVHC has an anti-static chamber to provide consistent medication availability. It can be used right out-of-package saving valuable time and resources.

A full range of products is available, designed specifically for each patient type.

The Monaghan Advantage

- ❖ The original valved holding chamber design
- ❖ Cost effective
- ❖ Anti-static, can be used directly out of package
- ❖ Optimized delivery system for consistent medication availability
- ❖ Easy to clean - dishwasher safe
- ❖ Low resistance valve
- ❖ Recommended by leading MDI manufacturers
- ❖ Compatible with most commonly prescribed MDIs

THE FACTS

The AEROCHAMBER PLUS® Z STAT® aVHC continues to offer the best value for optimal aerosol drug delivery and economic value for hospitals.

Different chambers deliver different amounts of medication.¹ The AEROCHAMBER PLUS® Z STAT® aVHC provides a consistent dose of medication substantially equivalent to the MDI dose.²

INCREASES
AEROSOL
SUSPENSION
TIME



518.561.7330 phone
800.833.9653 toll free
518.561.5088 fax

Monaghan Medical Corporation
Plattsburgh, New York 12901
monaghanmed.com

EZ FLOW Exhalation Valve minimizes re-breathing and directs exhaled flow away from the patient's face.

FLOWSIGnal® Whistle alerts patients to breathe slowly (on mouthpiece and large mask).



Low Resistance Inhalation Valve allows patients to easily inhale their medicine.

Anti-Static Chamber can be used directly out of package and provides consistent medication availability.

Universal backpiece fits most commonly prescribed metered dose inhalers.

THE FACTS (continued)

Designed to respond to shallow or low volume breathing patterns, the extra-sensitive inhalation valve makes it easier for patients to inhale the medication.

Can be used with commonly prescribed Metered Dose Inhalers

REFERENCES

1. Drug Safety Update – Inhaled products that contain corticosteroids, July 2008. Medicines and Healthcare Products Regulatory Agency.
2. TMI Lab Data. Studies on file.

CLEANING IS EASY

- AEROCHAMBER PLUS® Z STAT® aVHC (12 month, reusable single patient use)
 - Designed for 12 months of use at home, clinic, or hospital environment
 - AEROCHAMBER PLUS® Z STAT® aVHC chambers are dishwasher safe – follow care and use instructions included with the chamber provided by your healthcare provider.

ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
79710	AEROCHAMBER PLUS® Z STAT® aVHC with Mouthpiece	cs/10
79750	AEROCHAMBER PLUS® Z STAT® aVHC with Mouthpiece	cs/50
88710	AEROCHAMBER PLUS® Z STAT® aVHC with Small Mask	cs/10
78710	AEROCHAMBER PLUS® Z STAT® aVHC with Medium Mask	cs/10
80710	AEROCHAMBER PLUS® Z STAT® aVHC with Large Mask	cs/10

Contact your Monaghan Sales Territory Manager for more information or a demonstration.

